

Brampton National Platelet Immunology Reference Laboratory (NPIRL) Platelet / HLA Testing Services

Post Transfusion Purpura Investigation NPIRL PL-06

TEST DESCRIPTION

Post Transfusion Purpura (PTP) is a rare delayed transfusion reaction where a patient develops a sudden and dramatic drop in their platelet count typically between 6-10 days following red blood cell transfusion. Anti-HPA-1a antibody is most commonly associated with this reaction, however other HPA antibody specificities have also been found. Investigation for suspected PTP includes HPA antibody testing and HPA genotyping.

SPECIMEN AND REQUISITION REQUIREMENTS

Specimen(s)

- Three (3) full 5-7 ml EDTA (lavender) tubes mixed thoroughly by gentle agitation <u>and</u> one (1) 10mL SST serum tube inverted 5-10 times to activate clotting.
- Label specimen(s) with the required minimum information: patient's last name, first name, PHN/ (PHIN) or hospital name or other unique identifier, date of collection and facility name.

Complete Requisition (must include)

- Patient's last name, first name, date of birth and PHN/(PHIN) or hospital number or other unique identifier
- Facility Name
- · Physician/Health Care Provider name
- Phlebotomist name, classification, initial
- · Date/time of collection
- · Name, facility, address, fax and contact number of individual to whom the report will be sent

Requisition(s)

Platelet Immunology Requisition_NPIRL

PRE-SHIPPING STORAGE

- If unable to ship within 24 hours of collection store refrigerated at 1-10°C. Refrigerated samples must arrive within 48 hours of collection.
- If unable to ship within 48 hours of collection; separate and freeze <u>serum</u> sample only.
- · Do not open EDTA Samples

SHIPPING INSTRUCTIONS

• Contact the National Platelet Immunology Reference Laboratory (NPIRL) prior to shipping samples so that information may be obtained to assess the priority of the requested investigation.

Shipping

- Ship in a container with an ice pack (gel pack stored at 2 8°C may be used for refrigerated samples).
- Select shipping method for container to arrive at testing site within 24 hours.
- Fax shipping information (way bill #) to 905-494-8149. Indicate the patient name and type of investigation on the fax.



Brampton National Platelet Immunology Reference Laboratory (NPIRL) Platelet / HLA Testing Services

Post Transfusion Purpura Investigation NPIRL_PL-06

SEND TO

Canadian Blood Services Brampton Centre National Platelet Immunology Reference Laboratory (NPIRL) 100 Parkshore Drive Brampton, ON L6T 5M1 Tel: 905-494-5257

Fax: 905-494-8149