

OCTAPLASMA: BLOOD BANK COMMUNICATION

Administration of Octaplasma® must be based on ABO-blood group. Accordingly, Octaplasma is available in 4 specific blood groups, all of which can be ordered through Canadian Blood Services.

Please see below for additional information:

Is there a need to update the Information System of blood bank?

- Yes. Octaplasma has two approved manufacturing sites, so there may be a need to update your Information System to account for two additional Facility Identification Numbers (FIN), as per ISBT barcoding (see below table for Octapharma FIN codes)
- Depending on your handling procedures, each bag can be scanned into your system using the ISBT barcoding on the labeling.

ISBT-128 UNIQUE IDENTIFIER

Every bag of Octaplasma is individually and uniquely numbered, and is marked with an ISBT code. The ISBT 128 donation identification number (data structure 001) is used for this purpose.

Its format is:

=αpppppyynnnnnnff



OCTAPLASMA: BLOOD BANK COMMUNICATION

Structure	Description	Value	Remarks
=	Data identifier	-	Fix
αpppp	Facility identification number (FIN)	X0001	Indicates manufacturing facility. For Octaplasma you can expect the following: X0001 = manufactured at Octapharma Vienna X0002 = manufactured at Octapharma Stockholm X0003 = manufactured at Octapharma Stockholm
yy	Year of production	e.g. 07	
nnnnnn	Individual Bag Serial Number	000001 – 999999	Starting with 000001 at beginning of every year
ff	Flag	60 - 96	Checksum (not for use by end user)
☒		A – Z; 1 - 9 (1 digit only, alphanumeric values)	Checksum (not for use by end user)